



**Burhani Engineering Technology**  
**Pressure and Vacuum**

ENGINEERING SAFER TOMORROW

# Explosion Proof Acrylic Vacuum Glove Box

Safe, Transparent & Controlled  
Vacuum Workstation

Manufactured by

**Burhani Engineering Technology**



SAFER  
OPERATIONS



FULL  
VISIBILITY



RELIABLE  
CONTAINMENT

Model No. **BET-EXPGB-500**



High-Clarity  
Acrylic Enclosure



Explosion-Proof  
Vacuum Configuration



Side Transfer  
Chamber



Secure  
Quick-Lock Clamps



Excellent  
Product Visibility



Laboratory-Ready  
Construction

## TECHNICAL OVERVIEW

|                  |                                          |
|------------------|------------------------------------------|
| Model No.        | BET-EXPGB-500                            |
| Product Type     | Explosion Proof Acrylic Vacuum Glove Box |
| Material         | Acrylic                                  |
| Glove Ports      | 2 Ports                                  |
| Transfer Chamber | Integrated Side Antechamber              |
| Construction     | Transparent High-Clarity Acrylic Chamber |
| Sealing System   | Multi-point Quick-Lock Clamps            |
| Monitoring       | Pressure Gauge / Vacuum-Compatible Setup |
| Visibility       | Full Transparent Viewing Area            |

## APPLICATIONS

-  Chemical handling and processing
-  Solvent-sensitive work
-  R&D laboratories
-  Sample transfer
-  Controlled atmosphere operations

CONTAIN  
CONTROL  
EXPLORE  
SAFELY

A TRUSTED PARTNER  
IN VACUUM AND  
PRESSURE SOLUTIONS



### BURHANI ENGINEERING TECHNOLOGY

Plot no XVI, MAHAVEER ESTATE,  
CHETTIMEDU VILLAGE,  
PAYASAMABAKKAM, MADHAVARAM TALUK,  
CHENNAI - 600052



Tel: 9677296252



GSTIN: 33ACNPH2765G1ZM

Website QR



<https://bhfvacuum.com/>

Email QR



[mailto:sales@bhfvacuum.com?](mailto:sales@bhfvacuum.com?subject=Product%20Enquiry)  
subject=Product%20Enquiry

IDEAS  
ENGINEERED  
FOR A  
SAFER WORLD

VACUUM SOLUTIONS  
FOR A SAFER,  
CLEANER TOMORROW

*Precision  
Protection  
Possibilities*

