



Burhani Engineering Technology
Pressure and Vacuum

ENGINEERING
CONTROLLED
ATMOSPHERES
FOR A CLEANER TOMORROW

PRECISION
SAFETY
RELIABILITY
FOR YOUR RESEARCH

Nitrogen Purging Acrylic Glove Box

Model: BET-LAGB-A11

Precision-engineered acrylic glove box with integrated antechamber for controlled atmosphere handling and nitrogen purging applications.



High-transparency
acrylic chamber



Integrated
antechamber /
pass box



Dual glove ports
for controlled
handling




Pressure gauge
for monitoring



Gas inlet / outlet
provision



Controlled
atmosphere
operation



IDEAL FOR
R&D, RESEARCH LABS
AND CONTROLLED
ENVIRONMENT

SEE A CLEARER TOMORROW

INNOVATIVE
SOLUTIONS
FOR A CLEANER
SAFER WORLD

Dimensions & Spatial Specifications

Chamber Component	Metric Dimensions	Imperial Dimensions
Main Box Chamber (Internal)	700 W × 450 D × 550 H mm	27.5" W × 17.7" D × 21.6" H
Antechamber / Pass Box	240 W × 240 D × 240 H mm	9.4" W × 9.4" D × 9.4" H
Total System Footprint (Overall)	~1000 W × 750 D × 700 H mm	~39.4" W × 29.5" D × 27.5" H
Glove Port Diameter	Ø 190 mm	Ø 7.5"

CONTROL
CONTAIN
PURGE
PERFORM

BUILT FOR
A CLEANER,
SAFER TOMORROW



BURHANI ENGINEERING TECHNOLOGY

Plot no XVI, MAHAVEER ESTATE,
CHETTIMEDU VILLAGE,
PAYASAMBAKKAM, MADHAVARAM TALUK,
CHENNAI - 600052



Tel: 9677296252



GSTIN: 33ACNPH2765G1ZM



Website QR



<https://bhfvacuum.com/>



Email QR



[mailto:sales@bhfvacuum.com?
subject=Product%20Enquiry](mailto:sales@bhfvacuum.com?subject=Product%20Enquiry)