




Burhani Engineering Technology
Pressure and Vacuum

BET-NGB-S

SINGLE HEPA FILTRATION

Single HEPA filtration system provides HEPA filtered air into OR out of the glove box.



HEPA Filtered
Inlet Air



Negative
Pressure Chamber



Exhaust Through
HEPA Filter



Pressure
Monitoring




Safe. Ergonomic. Efficient.
Engineered for Maximum Protection.

BET-NGB-D

DOUBLE HEPA FILTRATION

Double HEPA filtration system provides HEPA filtered air into AND HEPA filtered air out of the glove box.



HEPA Filtered
Inlet Air



HEPA Filtered
Exhaust Air



Two (2) 8"
HEPA Filters



Enhanced
Contamination
Control



Dual Protection. Clean Air In & Out.
Superior Filtration for Critical Applications.

BET-NGB-CLC

CLOSED LOOP DOUBLE HEPA FILTRATION

Closed-Loop, Double HEPA filtration system continuously recirculates the internal atmosphere through two (2) 8" round HEPA filters.



Closed Loop
Recirculation



Dual HEPA
Filtration



Continuous
Internal Circulation



Maximum
Containment



Continuous Protection. Maximum Safety.
Closed-Loop Filtration for Highest Security.